



V300 30°

Lightweight wheelchair with reclinable backrest

The V300 30° is a lightweight manual wheelchair for temporary or permanent use with many adjustment possibilities. Changing positioning during the day, the use or the evolution of the mobility of the user, has many advantages. The V300 30° gives you the option of reclining the backrest.



VERMEIREN
we care for you

Quality guaranteed

Regarding quality, the **V300 30°** complies with the highest norms, laid down by certified test institutes in most European countries.

Colours

C59

champagne



Products designed to suit you

Our basic principle is: the wheelchair must be adapted to you, not the other way around. That is why our wheelchairs are evolutive and an answer to your specific needs.

All types of wheelchairs can be adapted, have numerous options and can be fully personalized so that there is a wheelchair for every person. This is the result of our continuous research and experience since 1957.



3 armrests heights
different armrests depths



Seat and backrest
tension adjustable

Technical data

■ Total width	600/630/650/670/690/710 mm
■ Total height	850 - 940 mm
■ Total length	1040 mm
■ Folded width	310 mm
■ Seat width	390/420/440/460/480/500 mm
■ Height armrests	220 - 240 mm

■ Seat height	440 - 530 mm
■ Height backrest	420 mm
■ Seat depth	440 - 520 mm
■ Weight	18,3 kg
■ Maximum user weight	130 kg
■ Additional info	CE

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

VERMEIREN Group N.V.

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



For more information, contact your local dealer

VERMEIREN
we care for you