

# Eclips XXL

## XXL WHEELCHAIR

Aluminum foldable wheelchair with double cross. The Eclips XXL provides a natural comfort with its anatomical seat and backrest. This wheelchair, equipped with reinforced rear wheels, can support a maximum weight of 200 kg.

Easy transfers thanks to the foldable aluminum footplate in one piece.



**C27**  
blue

**200KG**



B78 - Anti-tipping




B74 - Brakes



**ECLIPS XXL**

Wide range of accessories and options.



850	375	600	200	530	620	1070	530	1140	38,5	200
900	375	650	200	530	620	1070	540	1140	39	200
950	375	700	200	530	620	1070	550	1140	40	200
1000	375	750	200	530	620	1070	560	1140	40,5	200



Sizes in mm, weight in kg. Deviations in sizes and weight (15 mm, 1,5 kg and 1,5°) are possible.