



# Inovys II *Evo*

Tilt-in-space wheelchair





## Maximum comfort

When designing a tilt-in-space wheelchair, our priority is comfort. Based on our successful Inovys range, this wheelchair also provides excellent anti-decubitus qualities. The design of

Inovys II Evo is smartly thought to suit individual needs. Our many years of experience ensure versatility, endless adjustment possibilities, and comfort on top of a refreshed design.

## Customizable adjustments

Throughout the day, the use, and the user evolution, adjustment possibilities are a must fit needs, and bring comfort. Our the tool-free system allows for changes to provide comfort to both the user and the caregiver.



## Tilt in space

Pressure relief is an important benefit of Inovys II Evo. The tilt-in-space system rotates smoothly around the gravity center of the user without stretching the skin of the user. Other benefits include better blood circulation, less strain on joints or muscles, and better spinal health.



*Push handles adjustable in height without additional tools*



*Armrests removable and adjustable in height and depth  
Adjustable seat depth*



*Highly adjustable legrests  
Inclinable footplate  
Calf support adjustment*



## Designed with care

Inovys II Evo offers excellent comfort with its pre-shaped adjustable seat and backrest cushions. Other cushion shapes and configurations are also available to guarantee individual adjustment for each user.

## Your story, your configuration

Not only does the chair allow it to vary its current configuration and position. Inovys II Evo can be adapted to suit the users' needs in the best way possible with our wide range of options and accessories. To accommodate new needs, individual elements can be interchanged even after the first configuration. A chair created for and by you.



*Available as E-Version (electric adjustable seat and backrest)*



*T30 version - Available as a transit version with T30 rear wheels*

## Available color



**C89 - Copper**



*B74 - Compatible with brakes (and lock) for the attending person*



*B66MAL - Optional adjustable hemiplegic armrest for the left or the right side*

## Technical data

■ Total width	630 - 680 - 730 mm
■ Total height	990 - 1060 mm
■ Total length	1010 - 1050 mm
■ Seat width	400 - 450 - 500 mm
■ Seat inclination	0° - 21°
■ Backrest inclination	0° - 54°

■ Seat height	460 - 530 mm
■ Height backrest	550 - 585 mm
■ Seat depth	460 - 520 mm
■ Weight	36 kg
■ Maximum user weight	135 kg
■ Additional info	<b>CE</b>

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

## VERMEIREN Group N.V.

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



**For more information, contact your local dealer**

**VERMEIREN**  
we care for you