



# V-Drive

Electrical assistant for  
your manual wheelchair

**V-drive** is an easy-to-use assistant for the carer to help them overcome uneven ground or to climb slopes. The push assistant provides support when needed.



**VERMEIREN**  
we care for you

## Comfort

The V-Drive system simplifies your movements with a manual wheelchair. The V-Drive takes over all the work from your pushing assistant. Taking a steep hill or inclinations, bumps in the roads or just a small journey outside... all this is just ideal for the V-Drive. A userfriendly controller ensures optimal quality and comfort for daily usage of your wheelchair in combination with the V-Drive.

We are offering two versions of the V-Drive:

### VERSION STANDARD

max. user weight 135 kg  
(wheelchair included)

### VERSION H.D.

max. user weight 200 kg  
(wheelchair included)



The height adjustment allows the motor to be positioned horizontally.

## Technical data

■ Seat width	390 - 480 mm
■ Speed	6 km/h
■ Driving range	15 km
■ Battery / motor	2 x 12V 12 Ah - 150 W
■ Weight	22,4 kg**
■ Maximum user weight	135 kg*

■ Seat width	500 - 600 mm
■ Speed	6 km/h
■ Driving range	15 km
■ Battery / motor	2 x 12V 12 Ah - 200 W
■ Weight	23 kg**
■ Maximum user weight	200 kg*

Measurements in mm and kg. Deviations in sizes and weight (15 mm and 1,5 kg) are possible  
Subject to technical alterations. Images for illustration purposes only. Certain features may be optional. \* incl. wheelchair \*\* incl. batteries

## Installation

A very simple installation is available for left or rightside mounting of the V-Drive. Disassembly of your V-Drive for transport inside your car can be done in a minute. Disconnect the battery pack and the driving system and you are ready to go.

If you prefer not to use your V-Drive along the way, it is sufficient to lift up the driving system, hang it to on to one of the pushing handles and you have regained manual control over your wheelchair.

## Products designed to suit you

**VERMEIREN** always strives to come up with the best solution available to encompass your personal needs. As technology advances, so do we and as such we can offer you the most up to date products. This is the result of our continuous research and experience since 1957.

## Quality guaranteed

Regarding quality, the **V-Drive** complies with the highest norms, laid down by certified test institutes in most European countries.

## VERMEIREN Group

Vermeirenplein 1-15  
B-2920 Kalmthout  
Belgium

[www.vermeiren.com](http://www.vermeiren.com)  
[info@vermeiren.com](mailto:info@vermeiren.com)



For more information, contact your local dealer

**VERMEIREN**  
we care for you