

SCOOTERS

NEW EDITION



VERMEIREN
we care for you

- 01 eris
- 02 ceres
- 03 carpo
- 04 carpo_{xd}
- 05 PTIONS

CRAFTED IN EUROPE

ELEVATING QUALITY AND FLEXIBILITY

At **VERMEIREN**, we take immense pride in the meticulous craftsmanship of our Eris, Ceres and Carpo scooter line, produced right here in Europe. This isn't just a manufacturing choice; it's a testament to our commitment to excellence, quality, and reliability.

EXPERTISE AND EXPERIENCE

Each scooter in our lineup is crafted with unwavering attention to detail and quality. Our production process, carried out by our skilled workforce, is rooted in decades of expertise. From the initial design phase to the finishing touches, every step is meticulously executed to ensure that our scooters not only meet but exceed the highest standards.

LOCAL PRODUCTION, GLOBAL IMPACT

When you choose Eris, Ceres, Carpo scooter, you're not just purchasing a mode of transportation; you're investing in European craftsmanship, quality, and reliability. Join us in experiencing the difference that comes from a product crafted with passion, expertise, and a commitment to excellence.



eris

The Eris scooter is your in- and outdoor scooter for short distances to the shop, friends or family. Complete with indicators and lighting in LED version. The shopping tray can be taken off so you can carry it with you when purchasing small items at your local grocery store.

6 km/h 10 km/h



DISPLAY
Easy to use display




SHOPPING BASKET
Detachable



STEERING COLUMN
Adjustable for greater driving comfort



LED LIGHTING SYSTEM
Front and rear LED flashing light system with direction indicators



SUSPENSION
Integrated shock absorption for a smooth driving experience



PIVOTING SEAT
Equipped with a large pocket at the back

MAX. SPEED

Up to
10 km/h

DRIVING RANGE

Up to
37 km

USER WEIGHT

Up to
140 kg



TECHNICAL INFORMATION

Total width	570 mm
Total length	1185 mm
Total height	1150 mm
Seat height	450-500 mm
Turning radius	2900 mm
Driving range*	30-37 km
Motor	320 W
Speed*	6 or 10 km/h
Weight	75 kg
Max. user weight	140 kg
Batteries	2x12V - 38Ah/50Ah

* The driving range is calculated according to European standards.

* Speed and driving range may deviate (ca. 15%).

Permitted dimensional and weight tolerances: 15 mm, 1,5 kg and 1,5°. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

ceres

Ceres SE is an outdoor scooter that harmoniously blends design, quality, and comfort, offering the optimal balance for everyday trips in terms of usability, comfort, and performance.



SUSPENSION
Integrated shock absorption for a smooth driving experience



PIVOTING AND ADJUSTABLE SEAT
For increased comfort, equipped with headrest



DISPLAY
Easy to use display



SHOPPING BASKET
Detachable



STEERING COLUMN
Adjustable for greater comfort



LED LIGHTING SYSTEM
Front and rear LED flashing light system with direction indicators



MAX SPEED

Up to
15 km/h

DRIVING RANGE

Up to
48 km

USER WEIGHT

Up to
140 kg



TECHNICAL INFORMATION

Total width	620 mm
Total length	1300 mm
Total height	1280 mm
Seat height	450-500 mm
Turning radius	2800 mm
Driving range*	30-48 km
Motor	450 or 650 W
Speed*	10 or 15 km/h
Weight	106 kg
Max. user weight	140 kg
Batteries	2x12V - 38Ah/50Ah

* The driving range is calculated according to European standards.
* Speed and driving range may deviate (ca. 15%).

Permitted dimensional and weight tolerances: 15 mm, 1,5 kg and 1,5°. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

carpo

Carpo SE is a modern outdoor scooter designed for daily use, providing best comfort and performance for extended rides in any weather condition. It boasts a fully equipped, sporty appearance.



ADJUSTABLE ARMRESTS

Adjustable in height and width for increased comfort



PIVOTING SEAT

Adjustable for increased comfort, equipped with a removable headrest



LCD DISPLAY AND USB CONNECTION

USB connection for smartphone charger easy to use LCD display



SHOPPING BASKET

Detachable



HIGH GROUND CLEARANCE

For best outdoor driving experience



LIGHTING SYSTEM

Front and rear lights with indicator



MAX SPEED

Up to
15 km/h

DRIVING RANGE

Up to
45 km

USER WEIGHT

Up to
150 kg



TECHNICAL INFORMATION

Total width	660 mm
Total length	1535 mm
Total height	1390 mm
Seat height	480-555 mm
Turning radius	2000 mm
Driving range*	45 km
Motor	750 W
Speed*	10 or 15 km/h
Weight	144 kg
Max. user weight	150 kg
Batteries	2x12V - 80Ah

* The driving range is calculated according to European standards.

* Speed and driving range may deviate (ca. 15%).

Permitted dimensional and weight tolerances: 15 mm, 1.5 kg and 1.5%.
Subject to technical alterations. Images for illustration purposes only.
Certain features may be optional.

carpo^{xd}

The XD version is specially engineered to deliver the best performance and driving experience on mountain and hilly region roads. Our smartly programmed gear ratio provides exceptional torque and stability, ensuring you can confidently ride even during prolonged use in challenging conditions.

13
km/h



ADJUSTABLE ARMRESTS

Adjustable in height and width for increased comfort



PIVOTING SEAT

Adjustable for increased comfort, equipped with a removable headrest



LCD DISPLAY AND USB CONNECTION

USB connection for smartphone charger, easy to use LCD display



SHOPPING BASKET

Detachable



HIGH GROUND CLEARANCE

For the best outdoor driving experience



LIGHTING SYSTEM

Front and rear lights with indicator



MAX SPEED

Up to
13 km/h

DRIVING RANGE

Up to
57 km

USER WEIGHT

Up to
150 kg



TECHNICAL INFORMATION

Total width	660 mm
Total length	1535 mm
Total height	1390 mm
Seat height	480-555 mm
Turning radius	4020 mm
Driving range*	57 km
Motor	750 W
Speed*	13 km/h
Weight	144 kg
Max. user weight	150 kg
Batteries	2x12V - 100Ah

* The driving range is calculated according to European standards.

* Speed and driving range may deviate (ca. 15%).

Permitted dimensional and weight tolerances: 15 mm, 1.5 kg and 1.5°. Subject to technical alterations. Images for illustration purposes only. Certain features may be optional.

OPTIONS

VERMEIREN offers a wide range of accessories for your scooter. From crutch holder to scooter cover, optimise your scooter use with a selection of our options.

	Eris	Ceres SE	Carpo SE Carpo XD
• SW31	●	●	●
• B83		●	●
• SW06		●	●
• SW03		●	●
• SW83		●	●
• SH08	●	●	●
• SW09		●	●
• SW84		●	●
• MOVO		●	●



SW31
Crutch holder



B83
Oxygen Cylinder holder



SH08
Scooter cover



SW09
Windshield

S = STANDARD option. Subject to technical alterations.
Images for illustration purposes only.



SW06
Rollator holder



SW03
Topcase



SW84
Elevated footplate



SW83
Oxygen holder



mevo

evolution



THE ULTIMATE 4 SEASONS CANOPY FOR MOBILITY SCOOTERS



Developed in Canada

Developed and designed by Quebec engineers.

NO MATTER THE WEATHER

The Movo canopy protects you first of all from the sun, but also from all weather conditions. Combining protective fabric against the rain, windshield and sun visor, it also protects against cold and wind.

FOR YOUR SCOOTER

Thanks to an ingenious installation using the hitch under the seat, Movo canopy is compatible with most models of **VERMEIREN** mobility scooters, including Ceres and Carpo.

A CANOPY THAT TURNS INTO A CABIN IN AN INSTANT

Thanks to its retractable technology, it does not need to be lowered at all times. You can choose the position that suits you, depending on the weather. You therefore have the option: canopy and cabin option with the same product.



SUPER-FAST INSTALLATION

You can install or remove the Movo canopy in less than 10 seconds. It's ultra-quick.



EASY TO ASSEMBLE

A detailed assembly manual and an accurate assembly video are available on our website, so you will be able to assemble the canopy in less than 15 minutes. Once assembled, the installation can be completed in less than 10 seconds.



COMPACT DESIGN

The width of the Movo canopy is 74 cm, which allows it to fit through all openings, even small interior doors.



ULTRALIGHT

With a total weight of 8,7 kg (19 lbs), the Movo canopy can be installed easily by one person.



Contact your local supplier for more information.

VERMEIREN GROUP

Vermeirenplein 1-15
B-2920 Kalmthout
Belgium

www.vermeiren.com
info@vermeiren.com



VERMEIREN

we care for you