


Material			
No.	Plan	Article	Quantity
1	X005135	1797400	2
2	Y002782	1483354	1
3	Y020746	2004754	1
4	Y020128	1798200	1
5A	Y020130	1798154	1
5B	Y020179	1798554	1
6	Y020178	1798354	1
7	M8x20	1001580	1
8	M6x16 ISO7380	1008459	2
9	M10x1	1999935	2
10	M5x16 DIN7991	1907493	2
11	---	1005899	1

5A For Navix
5B For Express

E000914 / B

REMOVABLE CONTROLLER LE				Plan No.
				E000914 / B
Scale: 1:2	Drafted	Verified	Confirmed	
Date	15/01/2014	01/07/2014	01/07/2014	
Initials	SD	SD	PJ	