


INOVYS 2 EVO B23 SEAT FRAME				-
				Plan No.
Sheet 1/2				E00041297/B1.1
Scale:	1:4	Drafted	Verified	Confirmed
Date	12/07/2022	12/01/2024	12/01/2024	
Initials	DR	PM	AS	



Material			
No.	Plan	Article	Quantity
1	Nut M8 Black	2994009	2
2	X005575	2023600	2
3	Y011640	1655154	1
4	Y011652	1655254	1
5	Y011748	1659900	2
6	Y012123	1675800	2
7	YP06458	2734454 (35)	1
	Y011405	1657354 (40)	1
	Y011746	1663054 (45)	1
	Y011938	1663154 (50)	1
8	YP06455	2734354 (35)	1
	Y010813	1638754 (40)	1
	Y007691	1662854 (45)	1
	Y009834	1662954 (50)	1
9	M6x16 DIN912	2994065	8
10	M8 DIN126	1996214	1
11	M6 DIN125	2994414	12
12	Cap 20x20	1001612	2
13	M6 DIN985	2994063	1
14	Gas spring 400N	1005602	1
15	M6x35	2994131	1
16	M6x45 ISO7380	2994086	2
17	M8x40 ISO7380	2994309	2
18	Control screw M7	1900008	1
19	M8x60 ISO7380	2995491	4
20	L35 Black Soft 35	2995554	1
	L35 Black Soft 40	2995555	1
	L35 Black Soft 45	2995556	1
	L35 Black Soft 50	2995557	1
21	Board 35	2996475	1
	Board 40	1901287	1
	Board 45	1901288	1
	Board 50	1901289	1
22	M8x45	1915945	2

INOVYS 2 EVO B23 SEAT FRAME				-
				Plan No. E00041297/B1.1
Scale: 1:4	Drafted	Verified	Confirmed	
Date	12/07/2022	12/01/2024	12/01/2024	
Initials	DR	PM	AS	