


Detail A  
Scale 1:4



Head Side

Foot Side

Install the foot pedal (6) on the foot side of the bed

1905317

E000596 / G

Sheet 1/2

# BLOCKSYSTEM 970 ALOIS

Plan No. / Plannr.

E000596 / G

Scale / Schaal: 1:8	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd
Date / Datum	09/12/2009	04/01/2016	04/01/2016
Init. / Init.	JG	TT	PJ



Material / Materiaal			
No. / Nr.	Plan	Article / Artikel	Quantity / Aantal
1	Y015560	17229..	1
2	Y015561	17230..	1
3	Y015563	17231..	1
4	Y015564	17279..	1
5	Y015571	1731200	1
6	Y015584	17309..	1
7	Y016204	17308..	1
8	Y016328	1731374	1
9	X002897	1561774	1
10	---	1905247	4
11	M8 DIN934	1001671	1
12	M10 DIN985	1001683	4
13	M6x6 DIN916/45H	1998138	7
14	M8 DIN125	1003199	1
15	Ø12	1003783	6
16	M10x35	1009442	4
17	M6x06	1906051	12
18	M8x60 DIN444LB	1904677	1
19	X001016 (M8 DIN985)	1001682	2
20	M12 PA6.6	1993362	4
21	M8x25	1990501	2
22	Y016961	1739900	2

1905317

E000596 / G

Sheet 2/2

<b>BLOCKSYSTEM 970 ALOIS</b>				Plan No. / Plannr.
				E000596 / G
Scale / Schaal: <b>1:8</b>	Drafted / Getekend	Verified / Nagezien	Confirmed / Bevestigd	
Date / Datum	09/12/2009	04/01/2016	04/01/2016	
Init. / Init.	JG	TT	PJ	